

INSTRUCTIONS FOR USE

Cannula with Guide

CSK-6130, # CSK-6131

CAUTION: Federal Law (US) restricts this device to sale by or on the order of a physician.

PRODUCT DESCRIPTION

The Cannula with Guide is provided sterile, for single-use only. This product may be used in conjunction with the following devices provided separately:

(1) EPI-Sense® Coagulation System, sterile, single-use (under separate IFU)

The Cannula is manufactured PVC free and is not made with natural rubber latex.

PRODUCT FEATURES




Figure 1. Cannula (CSK-6131) with Guide Key Features

- [1] Torquer
- [2] Vacuum Port;
- [3] Cannula Cap
- [4] Guidewire (0.026" Diameter)

INDICATIONS

The Cannula with Guide is indicated for use in laparoscopic or general surgery to provide access for operative and diagnostic instrumentation in body cavity spaces.

CONTRAINDICATIONS

There are no known contraindications.

⚠ WARNINGS ⚠

Insertion or removal of the Cannula with Guide while guidewire is extended may cause potential patient harm. Always fully retract the guidewire into the Cannula with Guide lumen.

The cannula is provided sterile and is intended for single use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.

Inspect the device packaging prior to use. If any breach of the packaging is found, the sterility of the product cannot be assured, and the product should not be used to prevent patient infection.

Inspect the cannula and guidewire prior to use. Ensure cannula distal end and guidewire are smooth with no sharp edges. Sharp edge can cause potential patient harm. If sharp edge is found, device should not be used.

Care should be taken when inserting the Cannula with Guide. Applying excess force could cause potential patient harm. To reduce friction during insertion, lubricate the cannula with sterile saline.

Care should be taken when manipulating the cannula or guidewire. Always ensure no tissue is caught by the guidewire and brought into the Cannula with Guide lumen as this may cause altered hemodynamics or unintended tissue damage.

Care should be taken when manipulating the guidewire and any over-the-wire devices. Excessive forces may damage the cannula and/or guidewire.

Avoid excessive pulling on the torquer. Excessive pulling on torquer may damage the cannula and/or guidewire

⚠ CAUTIONS:

- Avoid over-rotating the cannula with Guide. Over-rotation can cause the vacuum tubing of the cannula to kink, reducing the cannula suction, thus causing reduced visibility.
- Avoid over-inserting the Cannula with Guide into patient body. Over insertion may reduce cannula suction.
- Inspect the device prior to use. If any damage is found, the function of the product cannot be assured, and the product should not be used.

POTENTIAL COMPLICATIONS OF THE COAGULATION PROCEDURE

- Blunt trauma to adjacent organs
- Infection
- Pericardial effusion
- Vessel injury
- Tissue perforation
- Hemodynamic instabilities
- Arrhythmias
- Thromboembolic complication
- Hernia
- Pneumothorax

REQUIRED EQUIPMENT AND SUPPLIES PROVIDED BY HOSPITAL

- Vacuum Tubing Set (Sterile)

RECOMMENDED EQUIPMENT

- 35 cm, 5 mm, and/or 10 mm Scope, depending on Cannula use
- 1000 mL 0.9% Normal Saline or Sterile Water
- EPI-Sense Coagulation Device – Refer to Instructions for Use for the EPI-Sense Coagulation Device when being used with the Cannula with Guide.

AT COMPLETION OF PROCEDURE

Remove cannula from tissue, disconnect all tubes, and discard cannula and tubing sets following local governing ordinances and recycling plans for disposal or recycling of device components.

TROUBLESHOOTING

Symptom	Action
Not able to draw vacuum	<ul style="list-style-type: none"> • Check vacuum connections. • Ensure vacuum regulator is set to -200 mmHg. • Examine cannula vacuum lumen for clots. • If clogged, flush vacuum lumen with sterile water. • If problem persists, replace cannula.
Unable to advance instruments within cannula lumen.	<ul style="list-style-type: none"> • Lubricate lumen with sterile water. • Ensure instruments are correct size.

Symbols

	Manufacturer		Catalog Number
	Consult Instructions for Use		Lot Number
	Sterile by irradiation		Use-By Date
	Do not re-use		Does not contain natural rubber latex
	Caution		Do Not Use if Package is Damaged
	Do Not Resterilize		Caution: Federal law (US) restricts this device to sale by or on the order of a physician or other licensed practitioner
	Model Number		Non-pyrogenic
 Transit temperature limit		 Transit humidity limit	

ABBREVIATIONS

IFU	Instructions for Use
-----	----------------------

LIMITED WARRANTY

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. AtriCure's sole obligation under this warranty is limited to the repair or replacement of this instrument. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure's control directly affect the instrument and the result obtained from its use. AtriCure assumes no liability with respect to instruments deliberately mis-used or those reused, reprocessed or re-sterilized and makes no warranties expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such mis-used or reused instruments. AtriCure shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the deliberate mis-use or re-use of this instrument.

DISCLAIMER

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

 AtriCure Inc.
7555 Innovation Way
Mason, Ohio 45040 USA
+1 866 349 2342
+1 513 755 4100

Page Intentionally Left Blank