**Ground Pad Adapter Cable Kit: CSK-2040**

**Instructions for Use**

AtriCure Incorporated  
7555 Innovation Way  
Mason, Ohio 45040 USA  
Customer Service:  
1-866-349-2342 (toll free)  
1-513-755-4100 (phone)

EC REP AtriCure Europe B.V.  
De entree 260  
1101 EE Amsterdam  
NL  
+31 20 7005560  
ear@atricure.com

For Use only with the nContact Coagulation System

**INDICATIONS:**

For use only with nContact Generator, nContact Coagulation Device and CSK-2000 RF Coagulation Cable (RF Cable).

The CSK-2040 Ground Pad Adapter Cable is provided non-sterile and should not be sterilized. See cleaning instructions section of this document.

**CONTENTS:**


**CSK-2040 Ground Pad Adapter Cable Kit contains:**

- 1 ea CSK-2040 Ground Pad Adapter Cable
- 1 ea Instructions for Use (IFU) Brochure

**INSTRUCTIONS:**

**CSK-2040**

1. Remove cable from package.
2. Before every use inspect cable for damage. Pay particular attention to the insulation of the cable. If any part of the cable is damaged, discard without using. See next page for schematic.
3. Connect ground pad cable to female end of ground pad adapter cable with “nContact Inc. CSK-2040” inscribed on connector.
4. Connect other end of CSK-2040 cable to the “Indifferent Electrode” receptacle on the front of the generator.

5. Follow instructions for use of generator and coagulation device.
6. To disconnect Ground Pad Adapter Cable from the generator, grasp by connector cover and pull back. Do NOT pull on wire.

7. To disconnect the Ground Pad cable from the Ground Pad Adapter cable, grasp the over-molded portion of the Ground Pad Adapter cable and the molded portion of the Ground Pad Cable and pull apart. Do NOT pull on the wires.
8. Excess pulling and flexing of the cable may damage the cable rendering it non-functional.
9. Store cable in safe dry area until next use.

**CLEANING INSTRUCTIONS:**

CAUTION: Ensure Isopropyl Alcohol (IPA) is completely dry before operating the unit.

CAUTION: Avoid caustic or abrasive cleaners

1. Disconnect the cable from the unit before cleaning.
2. If the cable is contaminated with blood and other body fluids, it shall be cleaned before the contamination can dry (within two hours of contamination).
3. The outer surfaces of the cable shall be cleaned with 70% - 90% Isopropyl Alcohol (IPA) wipes for a minimum of two minutes.
4. Pay attention to all area where fluids or soil may gather, such as under/around the connector or any tight crevices/grooves.
5. Dry the cable with a dry, white lint-free cloth.
6. Conduct a final confirmation of the cleaning process by visually inspecting the white cloth for remaining soil.
7. If soil remains on the white cloth, repeat steps 3 through 6.
8. After cleaning is complete, turn the unit on to perform Power On Self-Test (POST). If any errors are received, contact AtriCure to begin return process.

**DISPOSAL:**

Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

**WARNINGS :**

- For information concerning all warnings, precautions, and troubleshooting see nContact Coagulation System Radiofrequency (RF) Generator Unit Model CS-3000 Operators Manual.
- If necessary, clean the cable with 70%- 90% Isopropyl Alcohol (IPA) wipes to ensure that the cable does not cause an infection.

**CAUTION:** Cables to surgical electrodes should be positioned to prevent contact with patient or other leads.

**EXPLANATION OF SYMBOLS FROM PACKAGE LABEL:**

	Caution	REF	Model Number		Follow instructions for use
	Manufacturer		Expiration Date		Batch Code
	Instrument supplied Non-Sterile		CE mark, Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC.		
	Temperature Storage Limit		Humidity Storage Limit		
	Waste Electrical Electronic Equipment	<b>Rx ONLY</b>	Federal (US) law restricts this device to sale by or on the order of a physician.		

<table border="1"> <tr> <td data-bbox="44 95 142 159">EC</td> <td data-bbox="142 95 310 159">REP</td> </tr> </table>	EC	REP	Authorized Representative in the European Community
EC	REP		

RF=Radiofrequency                      LBL=Label


**WARRANTY INFORMATION:**

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure control directly affect the instrument and the result obtained from its use. AtriCure's obligation under this warranty is limited to the repair or replacement of this instrument and AtriCure shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AtriCure assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.


Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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