System Cart

User Manual

System-vogn - Brugervejledning (da)
Systemwagen – Benutzerhandbuch (de)
Carro del sistema – Manual del usuario (es)
Chariot du système – Manuel d'utilisation (fr)
Sistema carrello – Manuale d'uso (it)
Systeemkar – Gebruikershandleiding (nl)
Carrinho do Sistema – Manual do Usuário (pt)

P000880.M

AtriCure System Cart – User’s Manual
en

TABLE OF CONTENTS

1.0. Basic Cart	2
2.0. How to Use this Manual	3
3.0. System Description	3
4.0. Proper Use	4
5.0. EMC Guidance and Manufacturer’s Declaration	6

Rx ONLY

 **CAUTION:** Federal (U.S.) law restricts this device to sale by or on the order of a physician.

Manufactured by:
AtriCure Incorporated
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service:
1-866-349-2342 (toll free)
1-513-755-4100 (phone)

Information in this document is subject to change without notice. Reproduction of this document in any manner whatsoever without the written permission of AtriCure, Inc. is strictly forbidden. AtriCure and the AtriCure Logo are trademarks of AtriCure, Inc. AtriCure disclaims proprietary interest in the marks and names of others.

Printed in the USA.

1. SYSTEM OVERVIEW

1.0. Basic Cart

AtriCure developed a 'modular' approach to the cart. The System Cart consists of the Basic Cart (ASCB2) and two optional upgrade kits: (1) the Cryo Cart Kit (ASCC2); and (2) the ORLab Cart Kit (ASCO).

Figure 1: Basic Cart (ASCB2)

1.1. Customer Options

The modular approach allows customization to support all AtriCure components and accessories.


Figure 2: (1) Basic Cart; (2) Basic Cart with RF Generator; (3) Basic Cart with cryoICE BOX; (4) Basic Cart, RF Generator, Cryo Cart Kit and cryoICE BOX; (5) Basic Cart, RF Generator, ORLab Cart Kit and ORLab components; and (6) Basic Cart, RF Generator, ORLab Cart Kit, ORLab components, Cryo Cart Kit and cryoICE BOX.

2.0. How to Use this Manual

The following information is essential for safe and efficient operation of the AtriCure System Cart.

2.1. Table of symbols used in this manual

Item #	Symbol	Description
1		Consult operating instructions.
2		Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.
3		Manufactured by.
4		Alternating Current.
5		Temperature Storage Range
6		Humidity Storage Range
7		Voltage.
8		Waste Electrical and Electronic Equipment (WEEE)
9		Non-Sterile
10		Caution
11	Rx ONLY	CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

Table 1: Symbols

3.0. System Description

The System Cart supports medical equipment used in procedures for the treatment of arrhythmias in cardiothoracic theatres and hospital electrophysiological laboratories.

The complete system is comprised of five sub-systems:

1. **AtriCure® RF Generator System** produces and delivers RF energy, in a bipolar mode, at a frequency of approximately 460 kHz, with a maximum output power

ranging from 12 Watts up to 30 Watts depending on the operating mode. The system has an interface module that allows simultaneous connection of AtriCure® Isolator™ Transpolar™ ablative handpiece and pen devices to the RF generator.

2. **AtriCure® cryoICE BOX (ACM)** is an electro-mechanical and pneumatic cryogenic surgical system that delivers a cryogenic energy source, namely Nitrous Oxide, to create lines of ablation through cardiac tissue.
3. **AtriCure® ORLab** is a programmable diagnostic cardiac stimulator with electrocardiogram display tailored for use as a diagnostic tool during surgical and RF ablation procedures.
 - a. **NOTE:** Lexmark MS310D printer is used in conjunction with the ORLab to produce paper copies of diagnostic read outs.
 - b. **NOTE:** The handpiece and pen devices are separate components and are not a part of the AtriCure® Basic Cart.
4. **AtriCure® System Cart** is a 450 mm wide cart with an onboard Isolation Transformer. The Isolation Transformer supplies power to the other subsystems of AtriCure® System Cart.

No part of the AtriCure® System Cart or its accessories is intended to make contact the patient, instead, pacing stimuli are delivered, and ECG signal recorded, via AtriCure® Handpieces and approved third party products, which may include transvenous endocardial electrode catheters or epicardial mapping.

Indications for Use: The AtriCure System Cart is intended to be used as an accessory to the AtriCure Bipolar System, AtriCure Cryo Module System, and ORLab.

3.1. Equipment Classification

- 3.1.1. The AtriCure® System Cart consisting of the Basic Cart, the Cryo Cart Kit, and the ORLab Cart Kit are categorized as: Class 1 Equipment; Type CF Applied Part; IPX0; rated for 115-230 VAC 50-60 Hz; and 1800 VA.

4.0. Proper Use

- 4.1. **NOTE:** The AtriCure® System Cart consists of the Basic Cart and is ONLY intended to for use for the equipment in Table 2. A Cryo Cart Kit and an ORLab Cart Kit are also available and are ONLY intended for use with the equipment listed in Table 2.
- 4.2. **NOTE:** The cart should never be used to carry equipment other than what is specified in this manual. The cart should not be operated on an incline exceeding 10°. During operation and storage, the casters should be locked securely.

4.3. **NOTE:** This manual does not address specific instructions related to the equipment listed in Table 2. The user should refer to product specific manuals provided with individual components and systems for information pertaining to: (1) explanation of the function of controls, displays, and signals; (2) the sequence of operation; (3) the connection and disconnection of detachable parts and accessories; (4) indications of recognized accessories, detachable parts, and consumable materials; and (5) signal output or signal input ports intended only for connection to specified equipment described.

4.4. **NOTE:** Unplug cart before transport.

Manufacturer	Equipment Description	Product Code
MicroPace	ORLab EP Cardiac Stimulator/Recorder System	ORLab
AtriCure	AtriCure RF Generators	ASU2/ASU3/ASB3/MAG*
AtriCure	cryoICE BOX	ACM1/ACM2

*not all products are approved in all regions

Table 2: Compatible Equipment List.

4.5. Cleaning Instructions

4.5.1. **NOTE:** The Basic Cart Cleaning and Disinfection Instructions: (1) use a mild detergent (prepared to its specifications) and a damp cloth; and (2) do not use caustic, corrosive, or abrasive cleaning materials.

4.6. Serviceable Parts

4.6.1. **NOTE:** The System Cart consisting of: the Basic Cart; the Cryo Cart Kit; and ORLab Cart Kit does not have any serviceable parts. For servicing issues, contact AtriCure, Inc.

4.7. Environmental

4.7.1. **NOTE:** Environmental Conditions: (1) operational temperature +10° C to +40° C; (2) storage temperature -35° C to +54° C; and (3) humidity 15% to 90% relative humidity.

4.8. Preventative Maintenance

4.8.1. **NOTE:** Perform annual preventative maintenance procedures to ensure all Cart components are functioning properly, including but not limited to: (1) electrical power cords for fraying, damage, and proper grounding; (2) AC power switch; and (3) caster damage.

4.9. Safety Information

- 4.9.1. **⚠️ PRECAUTION:** To avoid electrical hazards, all parts of the ORLab, including the computer, monitor, Stimulus Generator Unit, and Lexmark laser printer must be connected to the Isolation Transformer and never directly to a mains outlet.
- 4.9.2. **⚠️ PRECAUTION:** To avoid risk of electric shock, the Isolation Transformer must only be connected to a supply mains with protective earth.
- 4.9.3. **⚠️ PRECAUTION:** To avoid electrical hazards and prevent overloading the transformer, do not connect any non-AtriCure approved equipment to the supplied Mains Isolation Transformer.
- 4.9.4. **⚠️ PRECAUTION:** Use AtriCure® System Cart only in ventilated areas and away from flammable gasses. To avoid risk of explosion, the Basic Cart should only be used in a ventilated area as gasses may be released during charging of backup battery, and should not be used in rooms with flammable anaesthesia.
- 4.9.5. **⚠️ PRECAUTION:** AtriCure® System Cart is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous.
- 4.9.6. **⚠️ PRECAUTION:** Secure tank with chain and ONLY use N2O Tanks ≤ 20 cm (8") in diameter that weigh ≤ 29.5 kg (65 pounds). Due to the weight of the N2O tanks, use caution while lifting and securing the tanks.
- 4.9.7. **⚠️ PRECAUTION:** If the cryoICE Box is connected to a grounded hospital AC power receptacle, it must be connected into the same grounded hospital AC power receptacle as Isolation Transformer.
- 4.9.8. **⚠️ PRECAUTION:** When the Isolation Transformer is plugged into the mains outlet, leave at least 2 ft clearance between the back of the AtriCure® System Cart and the wall to allow access for disconnection.
- 4.9.9. **⚠️ PRECAUTION:** AtriCure® System Cart atmospheric storage condition is 68.6 kPa to 101.3 kPa and the operating condition is 69.1 kPa to 101.3 kPa
- 4.9.10. **⚠️ PRECAUTION:** No modification of this equipment is allowed.

5.0. EMC Guidance and Manufacturer's Declaration

5.1. Electromagnetic Emissions

The AtriCure System Cart (ASC2), along with its associated equipment (ASU2/ASU3, ASB3, ACM1/ACM2, MAG and ORLab), from here on referred to as the AtriCure RF, Cryo, and EP Mapping System (ARC-EPMS) have been tested and found to comply with the limits for medical devices in IEC 60601-1-2:2014. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

The ARC-EPMS can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.

Guidance and manufacturer's declaration – electromagnetic emissions		
The AtriCure RF, Cryo, and EP Mapping System (ARC-EPMS) is intended for use in the electromagnetic environment specified below. The customer or the user of the ARC-EPMS unit should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ARC-EPMS unit, with the exception of the ASU or MAG RF Generator only during RF energy delivery, uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The ARC-EPMS unit is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

5.1.1. Electromagnetic Immunity

Guidance and manufacturer's declaration – electromagnetic immunity			
The AtriCure RF, Cryo, and EP Mapping System (ARC-EPMS) is intended for use in the electromagnetic environment specified below. The customer or the user of the ARC-EPMS unit should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Guidance and manufacturer's declaration – electromagnetic immunity			
Electrical fast transient/burst IEC 61000-4-4	± 2 kV @ 100 kHz repetition frequency for power supply lines ± 2 kV @ 100 kHz repetition frequency for input/output lines	± 2 kV @ 100 kHz repetition frequency for power supply lines ± 2 kV @ 100 kHz repetition frequency for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Power inputs ± 0.5 kV, ± 1 kV Line-to-Line ± 0.5 kV, ± 1 kV, ± 2 kV Line-to-Ground Signal input/outputs: ± 2 kV Line-to-Ground	Power inputs ± 0.5 kV, ± 1 kV Line-to-Line ± 0.5 kV, ± 1 kV, ± 2 kV Line-to-Ground Signal input/outputs: ± 2 kV Line-to-Ground	Mains power quality should be that of a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	0,15 MHz – 80 MHz 3V, 80 % AM at 1 kHz ISM bands between 0.15 MHz and 80 MHz 6V, 80 % AM at 1 kHz	0,15 MHz – 80 MHz 3V, 80 % AM at 1 kHz ISM bands between 0.15 MHz and 80 MHz 6V, 80 % AM at 1 kHz	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips: 0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° Voltage interruptions: 0 % UT; 250/300 cycle	Voltage Dips: 0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0 % UT; 1 cycle Single phase: at 0° 70 % UT; 25/30 cycles Single phase: at 0° Voltage interruptions: 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ARC-EPMS unit requires continued operation during power mains interruptions, it is recommended that the ARC-EPMS unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: <i>UT is the a.c. mains voltage prior to application of the test level.</i>			

5.1.2. EMC Guidance and Manufacturer's Declaration

Immunity test	Band (MHz)	Wireless Service	Immunity Test Level (V/m)	Compliance Test Level (V/m)
150 kHz to 80 MHz		General	< 3	< 3
	80 MHz – 2,7 GHz	General	3	3
	380 – 390	TETRA 400	27	27
	430 – 470	GMRS 460, FRS 460	28	28
	704 – 787	LTE Band 13, 1	9	9
	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	28	28
	1,700 – 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	28	28
	2,400 – 2,570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	28	28
	5,100 – 5,800	WLAN 802.11 a/n	9	9

Portable and mobile RF communications equipment should be used no closer to any part of the ASU or MAG System including cables, than the recommended separation distance calculated from the equation:

$$d = 6/E \times \sqrt{P}$$

Where:

d is the separation in meters

P is the maximum output power rating of the transmitter in watts (W) according to the Service

E is the Compliance Test Level indicated above.

Interference may occur in the vicinity of equipment marked with the following symbol:

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ASU System or any of its components are used exceeds the applicable RF compliance level above, the ASU System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating components or the entire ASU System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

5.1.3. Recommended Separation Distance

Recommended separation distances between portable and mobile RF communications equipment and the AtriCure RF, Cryo, and EP Mapping System

The AtriCure RF, Cryo, and EP Mapping System (ARC-EPMS) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ARC-EPMS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ARC-EPMS as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

5.1.4. Responsibility of the Manufacturer

AtriCure is responsible for safety, reliability, and performance of the equipment only if:

- Installation procedures in this manual are followed.
- Persons authorized by AtriCure carry out modifications or repairs.
- The electrical installation of the relevant room complies with local codes and regulatory requirements such as IEC and BSI.
- The equipment is used in accordance with the AtriCure User's Manual.

DISCLAIMER

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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