Sensing Cable: CSK-2030

Instructions for Use

AtriCure Incorporated
 7555 Innovation Way
 Mason, Ohio 45040 USA
 Customer Service:
 1-866-349-2342 (toll free)
 1-513-755-4100 (phone)

AtriCure Europe B.V.
 De Entree 260
 1101 EE Amsterdam
 The Netherlands
 +31 20 7005560
 ear@atricure.com

For Use only with the nContact Coagulation System

INDICATIONS:

For use only with nContact Generator, nContact Coagulation Device and CSK-2000 RF Coagulation Cable Extension (RF Cable).

The CSK-2030 Sensing Cable is provided non-sterile and should not be sterilized. See cleaning instructions section of this document.

CONTENTS:

CSK-2030 Sensing Cable Kit contains:

1 ea CS-2030 Sensing Cable 1 ea Instructions for Use (IFU) Brochure

INSTRUCTIONS:

CSK-2030

1. Remove cable from package.
2. Before every use inspect cable for damage. Pay particular attention to the insulation of the cable. If any part of the cable is damaged, discard without using. See next page for schematic.
3. Connect BLACK end of Sensing Cable "RF Generator" to BLACK port "Coagulation" on nContact Generator until a click is felt.
4. Connect BLACK end of CSK-2000 Cable to the BLACK receptacle "RF Cable" on CSK-2030 until a click is felt.
5. Connect 2mm shrouded pins to EP sensing equipment.
6. Follow instructions for use of generator and coagulation device.
7. To disconnect the RF cable from the CSK-2030 Enclosure, grasp the Enclosure with one hand, and with the other hand, grasp the connector cover and pull back. Do NOT pull on wires.
8. To disconnect Sensing Cable from the generator, grasp by connector cover and pull back. Do NOT pull on wire.
9. To disconnect 2mm shrouded pins, grasp the over-molded portion of the 2mm pin where the pin number is located and pull. Do NOT pull on the wire.
10. Excess pulling and flexing of the cable may damage the cable rendering it non-functional.
11. Store cable in safe dry area until next use.

CLEANING INSTRUCTIONS:

1. Disconnect the cable from the unit before cleaning.
2. If the cable is contaminated with blood and other body fluids, it shall be cleaned before the contamination can dry (within two hours of contamination).
3. The outer surfaces of the cable shall be cleaned with 70% - 90% Isopropyl Alcohol (IPA) wipes for a minimum of two minutes.
4. Pay attention to all area where fluids or soil may gather, such as under/around the connector or any tight crevices/grooves.
5. Dry the cable with a dry, white lint-free cloth.
6. Conduct a final confirmation of the cleaning process by visually inspecting the white cloth for remaining soil.
7. If soil remains on the white cloth, repeat steps 3 through 6.

DISPOSAL:

1. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

PRECAUTIONS:

1. Inappropriate handling of the connecting cable, including sterilization and immersing the electrical connectors, can result in degradation of system performance including the inability to initiate or complete and ablation therapy.

WARNINGS:

- For information concerning all warnings, precautions, and troubleshooting see nContact Coagulation System Radiofrequency (RF) Generator Unit Model CS-3000 Operators Manual.
- Ensure that the EP Sensing (EKG) equipment complies with IEC 60601-2-25 for protection from high frequency surgical interference.
- If necessary, clean the cable with 70% -90% Isopropyl Alcohol (IPA) wipes to ensure that the cable does not cause infection.

CAUTION: Cables to surgical electrodes should be positioned to prevent contact with patient or other leads.

EXPLANATION OF SYMBOLS FROM PACKAGE LABEL:

	Caution		Non-sterile		Lot Number		Storage Humidity Range
	Consult operator's manual		Waste Electrical Electronic Equipment		Catalog Number		Storage Temperature Range
	Manufacturer		Date of Manufacture		CE Mark and Notified Body		Authorized Representative in the European
Rx ONLY		Federal (US) law restricts this device to sale by or on the order of a physician.					

NOMENCLATURE:

RF=Radiofrequency LT=Label

WARRANTY AND DISCLAIMER INFORMATION:

AtriCure warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond AtriCure control directly affect the instrument and the result obtained from its use. AtriCure's obligation under this warranty is limited to the repair or replacement of this instrument and AtriCure shall not be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse of this product, including any loss, damage, or expense which is related to personal injury or damage to property. AtriCure neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AtriCure assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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