IMPLANT CARD LEAFLET

MD

CAUTION: Federal law (US) restricts this device to sale by or on the order of a physician.

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en

INSTRUCTIONS FOR COMPLETION:

The following information should be completed by the healthcare institution/provider on the front of the implant card. (See FIGURE 1).

?	Patient Identification
31	Procedure Date
+	Name and Address of the Implanting Healthcare Institution/Provider

Once the implant card is populated with the information above, peel the card from the backer, fold along the perforation, and bond it together. Give both the implant card leaflet and implant card to the patient.

PATIENT INFORMATION:

The implant card serves as a record of the LAA Exclusion System implant you have received. The information on the implant card helps communicate information about this implant. It can be presented to physicians or airport security. The implant card should be presented anytime you are subject to a Magnetic Resonance Imaging (MRI) screening, as the information helps ensure the MRI does not interfere with the implanted device.

In the event you need to get a replacement card for your implant, you will need to have the UDI and lot information for your device off the rear of the card (See FIGURE 2). Contact the manufacturer to request a replacement card.

If you do not have the implant card, contact your healthcare provider or consult your patient chart for UDI and lot information.

In the unlikely event of any issues with the implant, the information on the implant card will ensure that your physician can inform you if you are impacted.

Device lifetime is equal to patient lifetime.

MR CONDITIONAL

The LAA Exclusion System Clip is MR Conditional. A patient with the LAA Exclusion System Clip may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient:

- **Nominal Values of Static Magnetic Field:** 1.5-Tesla or 3.0-Tesla
- **Maximum Spatial Field Gradient:** 40-T/m (4,000-gauss/cm)
- **Type of RF Excitation:** Circularly Polarized (CP) (i.e., Quadrature Transmission)
- **Transmit RF Coil Information:** There are no transmit RF coil restrictions
- **Operating Mode of MR System:** Normal Operating Mode
- **Maximum Whole Body Averaged SAR:** 2-W/kg (Normal Operating Mode)
- **Limits on Scan Duration:** Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks)
- **MR Image Artifact:** The presence of this implant produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant is located in the area of interest

ATRICLEIP: ACH1, ACH2, PRO1, PRO2, LAA0	V-Clip: PROV, ACHV	ATRICLEIP MINI: ACHM, PROM
DEVICE CONTAINS		
ATRICLEIP	V CLIP	ATRICLEIP MINI
Titanium Grade 2	Titanium Grade 5	Titanium Ti64
Polyurethane	Polyethylene Terephthalate	Polyethylene Terephthalate
Nitinol	Titanium Dioxide	Titanium Dioxide
Polyethylene Terephthalate		Nitinol
Titanium Dioxide		

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