**CE**  
**0086**

**Instructions for Use**  
**Soft Tissue Atraumatic Retractor**  
**Catalog No's 400-210, 400-220, 400-230**

**STERILE**

## READ ALL INSTRUCTIONS BEFORE USE!

Federal (USA) law restricts this device to sale by or on the order of a physician.  
Single Use Device. Do not resterilize.

### PRODUCT DESCRIPTION

The ESTECH Soft Tissue Atraumatic Retractor (STAR) is a single-use sterile device. The STAR consists of a set of silicone straps with a customizable retraction system that allows the surgeon to configure the straps in a manner that provides the best access to the surgical site. The STAR is available in different sizes to accommodate different size incisions.

### INDICATION FOR USE

The ESTECH STAR is intended for use in cardiothoracic procedures to retract soft tissues adjacent to a sternotomy, thoracotomy, or other surgical port.

### WARNINGS AND PRECAUTIONS

- Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Do not use this device if package is damaged or opened.
- Do not use this device if it shows signs of damage, e.g., crimps, kinks, or crushed areas.
- This device is for single use only. **DO NOT REUSE. DO NOT RESTERILIZE.**
- Proper surgical procedures and techniques are necessarily the responsibility of the medical profession. These instructions are furnished for informational purposes only. Each surgeon must evaluate the appropriate use of this device, case by case, based on medical training, experience, and the type of surgical procedure employed.
- Dispose of used products in accordance with established hospital protocol for hazardous waste.

Verify all user connections of this device are secure prior to use. **USAGE INSTRUCTIONS**

**NOTE:** Use of the STAR requires that an incise drape (3M Ioban™ 2 or equivalent) be used to cover the thoracotomy site. The adhesive pads on the STAR may not adhere properly to bare skin.

1. Using aseptic technique, remove the STAR from sterile packaging and place in the sterile field.
2. Fold the STAR in half and insert into thoracotomy incision with silicone straps perpendicular to the long axis of the incision. Holding the ring folded with one hand, push the exposed portion of the ring with opposite hand until the ring completely inside the chest

**NOTE:** STAR ring diameter should be greater than incision length to minimize potential for ring to withdraw from incision under tension.

3. Once inserted into the incision, grasp the straps and push the device in fully to allow the ring to deploy, resuming its circular shape. Gently pull the straps outward and align to create a uniform profile, covering all soft tissue.
4. Peel the backing off of one of the adhesive pads and fasten to the chest. Peel the backing off of the strap opposite to the first strap installed, and fasten to the chest wall. Continue fastening the straps the chest in this manner. Opposite straps will have the same logo printed on them (either ESTECH or STAR logo).

**NOTE:** Ensure that the draping is dry in the area that the adhesive pad will be fastened to. Avoid contacting adhesive pad with wet materials, as this will reduce adhesion.

**WARNING:** Avoid excessive tensioning of the straps and/or sequential fastening of adjacent straps. This may cause ring migration or dislodging of the device from the incision.

5. When all straps have been fastened, verify that all adhesive backings which were removed have been collected and discarded appropriately.
6. ESTECH recommends the use of the SI retractor in conjunction with the STAR to retract the ribs and enlarge the surgical window.

**WARNING:** If the SI retractor is used, do not attempt to spread the incision larger than the diameter of the metal ring. Doing so may cause damage to the SI retractor blades.

**NOTE:** If the SI retractor is used, select blades that are appropriate for the chest wall thickness. If blades that are much longer than the chest wall thickness are used, the SI retractor may migrate when expanded.

7. When finished, cut each of the straps distal to the adhesive pad and pull one of the straps away from the chest and toward the center of the incision to release the ring and remove retractor from incision.

**NOTE:** Removal of adhesive pads from sterile drape will cause drape to tear.

## WARRANTY

AtriCure provides no warranties, either express or implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. AtriCure shall not be liable to any person or entity for medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind AtriCure to any representation or warranty with respect to the product.


The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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